



Study Summary



Influence of frequency and amplitude on the mucus viscoelasticity of the novel mechano-acoustic Frequencer™

June 13th , 2019

Author: Dalma Schieppati, Rémi Germon, Federico Galli, Marco Giulio Rigamonti, Marta Stucchi, and Daria Camilla Boffito

www.dymedso.com

Purpose of the study


To study a mechano-acoustic treatment device to help patients expectorate excess mucus in a non-invasive and enervating way. It was important to understand how the mucus rheology and hydration changes depending on the frequency and amplitude of the acoustic waves delivered to the sample using the Frequencer™.

Background

Cystic Fibrosis (CF) is an inherited genetic disease affecting 1/3200 Caucasians that is characterized by the buildup of thick, sticky mucus that can damage many organs, especially the lungs [1]. CF patients live with mucus in their lungs that is dehydrated and viscous, making it harder to expectorate. This is namely caused by an absence or malfunction of the gene that encodes for the CFTR protein, responsible for keeping the mucus hydrated by maintaining the water and ion homeostasis across the respiratory epithelia [2]. Mucus is composed of water, protein, carbohydrates, lipids and antimicrobial factors [3]. Mucins are the main constituents of mucus and are responsible for its rheology [4]. Healthy humans have a mucin content of 1% by weight where individuals with CF have a mucin content of about 4% by weight [5].



Methods



Synthetic mucin solutions were prepared in vitro using porcine stomach mucin to recreate the 1% by weight healthy condition (S1) and the 4% by weight CF patient condition (S2). 48 experiments were conducted by adding mucus to a plastic cylindrical reactor with an external diameter corresponding to that of the opening of the Frequencer™ adaptor. Four conditions were tested: mucus on its own, mucus combined with water, mucus combined with NaCl (0.5 g L⁻¹) and mucus combined with brine. Multiple parameters on the Frequencer were set in order to determine which condition, frequency (20 Hz, 40 Hz, 60 Hz) and intensity (50% and 100% amplitude) would lead to optimal mucus rehydration and rheology.

Results


The Frequencer™ proved to be effective in the homogenization of synthetic mucin solution in vitro in 20 minutes. A working frequency of 40 Hz and a 0.5 g L⁻¹ NaCl solution are the optimal operative parameters to obtain partial rehydration of mucus, regardless of the intensity selected, this is consistent with what is reported by patients using the device.

Conclusion

The optimal frequency for mucus rehydration was found to be 40 Hz, regardless of the intensity selected and of NaCl concentration.

$$WA = \mu + \beta_1 F + \beta_2 S + \beta_3 M \cdot F + \beta_4 M \cdot S + \beta_5 F \cdot S + \beta_6 M \cdot F \cdot S$$

Mass fraction (M) Frequency (F) Salt concentration (S)



References

1. A. Stephenson, et al., Survival comparison of patients with cystic fibrosis in Canada and the United States: a population-based cohort study, *Ann. Intern. Med.* 166 (2017) 537–546, <https://doi.org/10.7326/M16-0858>.
2. I. McIntosh, G.R. Cutting, Cystic fibrosis transmembrane conductance regulator and the etiology and pathogenesis of cystic fibrosis, *FASEB J.* 6 (1992) 2775–2782, <https://doi.org/10.1096/fasebj.6.10.1378801>.
3. M.P. Rogan, et al., Antimicrobial proteins and polypeptides in pulmonary innate defence, *Respir. Res.* 7 (2006) 29, <https://doi.org/10.1186/1465-9921-7-29>.
4. D.J. Thornton, K. Rousseau, M.A. McGuckin, Structure and function of the polymeric mucins in airways mucus, *Annu. Rev. Physiol.* 70 (2008) 459–486, <https://doi.org/10.1146/annurev.physiol.70.113006.100702>.
5. S. Kirkham, J.K. Sheehan, D. Knight, P.S. Richardson, D.J. Thornton, Heterogeneity of airway mucus: variations in the amounts and glycoforms of the major oligomeric mucins muc5ac and muc5b, *Biochem. J.* 361 (2002) 537–546, <https://doi.org/10.1042/bj3610537>.